

Jasco™ diaphragm vacuum pump

Model No. RS-68

Product Description

“Jasco” RS-68 diaphragm vacuum pump has features of more advanced design, higher efficiency, longer operating life. It is an all new product engaging in the new and high technology. It is mainly used in medicinal products analysis, tenuously chemical engineering, biochemical pharmacy, food examination, investigating and solving criminal case, etc. It is an ideal product used with the high-accuracy chromatogram apparatus, and it is the essential equipment in the laboratory.

Pumping speed: (L/Min)	20	Pump Head	1
Ultimate Pressure	250mbar	Voltage Rating	230Vac, 50Hz
Dimensions (L x B x H) (mm)	235×140×210	Temp of the body (°C)	<55
Power (w)	200	Weight (Kg)	7.5
Connection: Inlet (mm) Outlet (mm)	$\phi 6$ Silencer	Material: Diaphragm Valves Pump head	HNBR HNBR PTFE coated
Working Temp (°C)	7—40	Noise Level (dB)	<50



Powered by DIYTrade.com

Powered by DIYTrade.com

International Co. for Lab. Instruments

45 Shalhoub st. Ain Shams, Cairo, Egypt+ Post Code 11311



+202 24916139 or +202 24916139



+202 24953441