

# CIRCULATOR OVEN

ITEM	CIRCULATOR OVEN					
MODEL	MAO-42C	MAO-43C	MAO-62C	MAO-63C	MAO-45C	MAO-65C
Chamber Size DxHxW cm	40x40x45	40x40x45	50x60x50	50x60x50	40x40x45	50x60x50
Overall Size DxHxW cm	56x90x56	56x90x56	66x112x60	66x112x60	60x94x60	70x116x64
Capacity	72 L	72 L	150 L	150 L	72 L	150 L
Temperature Range	Room Temp. ~200°C	Room Temp. ~300°C	Room Temp. ~200°C	Room Temp. ~300°C	Room Temp. ~500°C	Room Temp. ~500°C
Material	Outer: Galvanized Steel Plate with High Temperature Finish. Inner: Stainless Steel Plate					
Circulator Fan	Whole new made high-temperature resist long shaft motor and louvered fan. Use for CED convection to exhaust steam, waste gas and volatile gas, also dry rapidly.					
Temp. Controller	Digital LED display      Precision: 0.5°C      Spread: ±2°C					
Control System	SSR non-contact control system, quiet & stable					
Material of Heat	Adopt imported mineral wool. No pollution.			Imported heat-resistance ceramic wool & mineral wool, double temperature maintenance.		
Adjustable PALINGS	2 palings (steel plates)		3 palings (steel plates)		2palings (steel plates)	3 palings (steel plates)
Power Supply	220v / 10A	220v / 15A	220v / 1 5A	220v / 18A	220v / 24A	220v / 27A
Optional Accessories	<ol style="list-style-type: none"> <li>PID Micro computer temp. controller</li> <li>Digital time keeper, start to clock when reach the setting temperature. When time's up, heating stops and alarm rings. 1 timer, 1 double timer, 1 buzzer, button switch, blow holes, control lines.</li> <li>West Germany supper high-temperature protector.</li> <li>Transparent window, above 5 mm double reinforced glasses.</li> </ol>					

